



BioTeke

BioTeke Corporation (Wuxi) Co.,Ltd



Saliva Pregnancy Test Kit HIGH ACCURACY

www.bioteke.cn

Intended use



The Saliva Pregnancy Test Kit is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin (HCG) in saliva to aid in the early detection of pregnancy.

Advantage

CLEANER AND MORE CONVENIENT


Saliva sampling, a novel way of testing.

No urine sampling, no "dirty hands".




TESTING AT ANY TIME

Convenient testing, not restricted by time or environment. No need to go to the toilets, just sit at the table and test.





HIGH SENSITIVITY

Can detect pregnancy earlier. Higher sensitivity of 5 mlu/mL.



SUPERIOR PERFORMANCE

Stable performance, with accurate and reliable results.



PORTABLE

Easy to carry, reagent card light and easy to save.




SIMPLE TESTING PROCEDURE

Easy to test, using saliva to detect pregnancy. Easy to understand, clear and easy-to-read test results, providing a comprehensive overview.




Components




Components	Loading quantity(Specification)	
	1 Test/Kit	5 Test/Kits
Test Card	1pc	5pc
Diluent Tube	1pc	5pc
Tube Cap	1pc	5pc
Saliva Collector	1pc	5pc

Test procedure



Results



BioTeke Corporation

BioTeke Corporation (Wuxi) Co., Ltd.

For more information, you can call BioTeke national service hotline:

(0510) 8332 3992

Or visit the official website of BioTeke:

www.bioteke.cn

Version: 1.0
Revision date: Mar 18, 2024